Norsk

Wear’N’Go

Anklestøtte

Funksjon:

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English
Wear’N’Go – Ankle Support
Fuction:

Provides light support, compression, warmth, and stability to the ankle. The material increases blood circulation and relieves discomfort. The product is well-suited for use in both sports and rehabilitation after injuries or ailments.

Instructions for Use:

- Hold the product and pull it up over the ankle.
- Pull the product so that the silicone pad surrounds the ankle, the toes remain free, and the circle on the product is positioned over the heel.
- To remove the product, grip the lower edge above the toes and pull it off the foot/ankle.

Sizes:

- Measure as described on the packaging.
- Determine the correct size using the size guide on the packaging.

Materials:
See the packaging.

Washing Instructions:

Follow the instructions on the packaging.

Important Information:

- Use the product carefully and appropriately.
- Do not modify the product inappropriately / The function of the product must not be altered or weakened.
- Do not use if you are hypersensitive to any of the product’s components.
- The product is not flame-resistant.
- Avoid temperatures above 120 °C / 248 °F.
- Use the product only for the purpose described above.
- Recommended usage time is a maximum of 3–4 hours at a time.
- Try the product to ensure it fits correctly.
- For questions regarding use, contact the manufacturer or your local supplier.
- All users are advised to consult healthcare professionals before use. Athletes should consult a coach, physiotherapist, or doctor before use.
- Discontinue use and contact a healthcare professional if pain, skin irritation, injury, or unusual reactions occur.
- Do not use the product directly on open wounds.
- The product is intended for individuals over 18 years of age.
- Use is at the user’s own risk.

Disposal:

Dispose of the product and packaging in accordance with applicable local or national environmental regulations.

Liability Disclaimer:

The manufacturer assumes no responsibility in the following cases:

- The product is not maintained as described in the instructions.
- The product is combined with components from other manufacturers.
- The product is used outside of the recommended conditions, areas of use, or environment.
- Injuries resulting from misuse of the product.

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