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English Wear'N'Go Arm Compression Sleeves Function:

Provides firm and even compression to the arm, helping to relieve discomfort, improve blood circulation, support muscle recovery, and reduce the risk of future injuries. The product is well-suited for use during sports activities and for recovery after training.

Instructions for Use:

- Grasp the product and gently pull it over your hand and up your arm.
- Adjust the sleeve so it sits evenly, covering the area from the wrist to the upper arm. Ensure the fabric lies flat without wrinkles.
- To remove, grip the top of the sleeve at the upper arm and carefully pull it down and off.

Sizing:

- Measure as described on the packaging.
- Find the correct size using the size guide on the packaging.

Materials:

See the packaging.

Washing Instructions:

Follow the instructions on the packaging.

Important Information:

- Use the product carefully and appropriately.
- Do not modify or alter the product's function.
- Not to be used by individuals sensitive to any of the product's components.
- The product is not flame resistant.
- Avoid temperatures above 120°C / 248°F.

- The product should only be used as described above.
- Recommended usage time is a maximum of 3–4 hours at a time.
- Try the product to ensure proper fit.
- For questions about use, contact the manufacturer or local supplier.
- All users are advised to consult a healthcare professional before use. Athletes should consult a coach, physiotherapist, or doctor.
- Discontinue use and consult a healthcare professional if pain, skin irritation, injury, or unusual reactions occur.
- Do not use directly on open wounds.
- The product is intended for individuals over 18 years of age.
- Use at your own risk.

Product Disposal:

Dispose of the product and packaging in accordance with local or national environmental regulations.

Liability Disclaimer:

The manufacturer is not liable in the following cases:

- The product is not maintained as described in the instructions.
- The product is combined with components from other manufacturers.
- The product is used outside recommended conditions, use cases, or environments.
- Injuries caused by improper use of the product.

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