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Engelsk
Wear’N’Go
Knee Support
Function:
 Provides support, compression, and stability to the knee. The material enhances blood circulation and alleviates discomfort. The product is well-suited for use in both sports and rehabilitation after injuries or ailments.

Instructions for Use:

- Hold the product as shown in figures 1 and 2.
- Pull the product over the knee so that the silicone pad surrounds the kneecap, as shown in figures 3 and 4.
- To remove the product, grasp the bottom of the metal springs on both sides and pull it off the knee.

Sizes:

- Measure as described on the packaging.
- Determine the correct size using the size guide on the packaging.

Materials:
 See the packaging

Washing Instructions:
 Follow the instructions on the packaging.

Important Information:

- Use the product carefully and appropriately.
- Do not modify the product improperly; its function must not be altered or weakened.
- Do not use if you are sensitive to any of the product's components.
- The product is not flame-resistant.
- Avoid temperatures above 120°C / 248°F.
- Use the product only for its intended purpose

as described above.

- Do not use the product for more than 3–4 hours at a time.
- Ensure the product fits properly before use.
- If you have questions about using the product, contact the manufacturer or your local supplier.
- It is recommended to consult healthcare professionals before using the product. Athletes should consult with a coach, physiotherapist, or doctor before use.
- Discontinue use and consult healthcare professionals if you experience pain, skin irritation, injury, or unusual reactions while using the product.
- Never use the product directly over open wounds.
- The product is intended for individuals over 18 years of age.
- Use of the product is at your own risk.

Disposal:

Dispose of the product and packaging in accordance with applicable local or national environmental regulations.

Liability:

The manufacturer assumes no responsibility in the following cases:

- Products not maintained as described in the instructions.
- Products combined with components from other manufacturers.
- Products used outside recommended conditions, applications, or environments.
- Damage resulting from misuse of the product.

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